

# Operative | Experience

Emergency Obstetrics Surgery (EOS™)

Procedural Training

C-celia™ Post Partum Hemorrhage

Training Simulator

**Care, Maintenance and Operating Instructions Manual**

[www.operativeexperience.com](http://www.operativeexperience.com)

410-287-5748


## Mobius Retractor

The use of a Mobius Elastic Retractor, while not required or provided, will extend the life of the model.



### Parts List

- |                            |                            |                   |
|----------------------------|----------------------------|-------------------|
| A. Back plate              | B. Quick release pin       | C. Fill tubes     |
| D. Squeeze pumps           | E. Extra screws            | F. Simulator bags |
| G. Fill bags (blood, lube) | H. Pressure sleeve         | I. Intestines     |
| J. Uterine/pelvic lock     | K. Screws                  | L. Uterus         |
| M. Abdominal insert        | Not pictured – screwdriver |                   |



(Post Partum Hemorrhage Simulator)

2

## Installing Uterus

**Lubricate uterus for insertion**




**Identify hole in abdominal lining**




**Insert blood tube from uterus  
into hole**



**Iteratively and gently pull tube through  
as you push uterus into position. Alternate  
from one side to the other. Lubricating tube  
will ease insertion.**



**Uterus shown in good position**



(Post Partum Hemorrhage Simulator)

## Installing Uterus (cont.)


Continue sliding uterus down while gently pulling tubes...



...until uterus is in place




Push uterus into pelvis




seat the uterus



Check that tube is not kinked



Snap connect uterine blood flow tubes



(Post Partum Hemorrhage Simulator)

## Installing Uterus (cont.)

Image of connected blood flow tubes




Gently pull tubes from back until taut




Tuck tubes behind tissue



Position uterine/pelvic anchor



Align anchor posts with holes



Attach with T9 torx screws. Do not cross threads or overtighten.



(Post Partum Hemorrhage Simulator)

## Adding Intestines and Abdominal Wall

Place intestines in abdominal cavity




Locate attachment points on abdominal wall ring



Location of placement of attachment points



Position abdominal wall, slide attachment points under location of placement.



Insert quick release pin through frame



To remove abdominal wall, push the blue button on the quick release pin, extract pin and slide abdominal wall down.

(Post Partum Hemorrhage Simulator)


6

## Infusing with blood or lubricant

Open back plate and locate fill ports




Connect fill tube with IV fill bag



Snap connect to fill port



Place IV fill bag into pressure sleeve



To infuse:

- Close (turn to right) flow control knobs
- Put internal reservoir hand-squeeze pump in deflate (left) position
- Pump pressure cuff until blood/lube IV bag is empty
- Detach from fill port
- Release pressure on cuff


(Post Partum Hemorrhage Simulator)

7




## Pressurizing Internal Fluid Reservoir

Turn flow control OFF (turn right)  
Set hand squeeze pump to “inflate”



Pump hand squeeze pump...



...until pressure gauge shows green stripe



(Post Partum Hemorrhage Simulator)

## Lubricating and Bleeding

To lube, open valve and distribute lube in abdominal cavity around uterus with hand through abdominal opening. Close valve.



For arterial blood flow, turn blood flow valve to the left. Note: a  $\frac{1}{4}$  turn will provide strong blood flow. There is a separate valve on each side of the simulator.



## Disassembly, Cleaning and Storage

### Removing fluids from internal reservoirs.

- Turn off flow control valves
- Pressurize internal reservoir
- Insert fill tube IV spike into empty IV bag
- Attach fill tube connector to fluid port to extract fluid

### Flushing internal fluid reservoirs

Follow the same process for infusing model with blood or lube, substituting water or a mix of water and mild soap

### Cleaning

All parts of the model can be cleaned with water or water with a mild soap. Do not store the model overnight without washing, and wash and dry thoroughly before storage.

Note the entire model can be hosed down.

### Storage

All fluid should be removed from the internal reservoirs, and the system should be flushed with water before storage. Surface areas, including the inside compartment housing the infusion ports, should be cleaned and dried before storage.

**Technical Support:**

\* Simulators must be serviced by qualified OEI technicians only\*

OEI Office:	410-287-5748 x 302
Chief Engineer:	443-373-6815
Chief Operating Officer:	443-691-5179
Chief Executive Officer, Robert Buckman, Jr., M.D.	484-557-8876

[www.operativeexperience.com](http://www.operativeexperience.com)

**\* This Manual includes confidential, unpublished property of Operative Experience Inc. Photocopying and distribution are strictly prohibited. All Rights Reserved 2014.**

**Disclosure**

**To prevent electronic and mechanical system failures, please comply with the instructions included in this manual.**

**Infuse manikin with supplied or approved simulated blood and fluids only.**

**\* FAILURE TO COMPLY COULD RESULT IN CATASTROPHIC SYSTEM FAILURE\***



## Post Partum Hemorrhage Trainer

Mobius Retractor Information.....	1
Parts List.....	2
Installing Uterus.....	3-5
Adding Intestines and Abdominal Wall.....	6
Infusing with Blood and Lubricant.....	7
Pressurizing Internal Fluid Reservoir.....	8
Lubricating and Bleeding.....	9
Disassembly, Cleaning and Storage.....	10
Technical Support and Disclosures.....	11