*...going one step further*

3B  
LASER NEEDLE

---

# 3B LASER NEEDLE

Gebrauchsanleitung  
Operating Instructions  
Mode emploi  
Istruzioni per l'uso  
Manual de instrucciones  
Instrução de funcionamento

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E Dear User,  
 N We are glad that you have selected an ultra-modern, quality medical product in the 3B LASER NEEDLE made by 3B Scientific,  
 G and would like to thank you for placing your trust in us.  
 L In order to achieve the success you are looking for and to ensure that you have many years of use from your  
 I 3B LASER NEEDLE, we would ask you to read the following usage information carefully.  
 S If you have any questions, please feel free to contact us on +49(0)40 73966222 or send an e-mail to  
 H stefan.baudis@3bscientific.com.  
 Specific therapy advice can be obtained from specialist publications or by attending user courses and workshops  
 (see also www.3b-laser.de).

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## Operating instructions

1. Scope of delivery	16	7. Maintenance and care	24
2. Installation	17	8. Safety information	24
3. Use	18	9. Accessories / Label and warning signs	25
4. Contraindications	18	10. Status, Warning and Fault Messages	26
5. Operation	19	11. Technical Data	28
6. Usage	22	Annex 1	86

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





W14204 1 Pair of patient goggles (Softcap)

## 2. Installation

- Attach the table clamp (5) to the table or the equipment trolley and make sure that it is secure. If you are fixing it to an equipment trolley, make sure that the trolley itself is stable enough. A wall bracket is available for fixing to a wall, if required.
- Insert the holding arm (2) into the holding device as shown in Fig. 1 and locate it in the required position with the butterfly screws. Tighten the butterfly screws, making sure that they are not over-tightened.
- Fix the Laser Unit (3) to the metal plate using the supplied round nuts, as shown in Fig. 2. Please note: Tighten in a counter-clockwise direction.
- Connect the holding arm cable to the Laser Unit (3), as shown in Fig. 2.

### Connecting the fiber optic cables (FOC) (8):

The FOC (8) are numbered. Remove FOC 1 from the packaging and screw it in to position 1 on the Laser Unit, as shown in Fig. 3. Make sure that the FOC is not twisted or knotted. **(Note the minimum bend radius of 1.5 cm.)**


Feed the FOC upwards through the corresponding hole (Fig. 4) and attach it to the front of the acrylic plate (Fig. 5). Then remove the remaining FOC from the packaging in numerical order and connect them as detailed above, until all 12 FOC are in place (Fig. 6).

**FOC should not be kinked or placed under tension. Make sure that the FOC can hang freely and that there is sufficient space to the floor. Do not tread on the FOC or apply any strong pressure to them. Only handle FOC using the white holders or the screw threads.**

- Connect the power supply unit (10) with the holding arm cable and plug it into the mains power (100 VAC - 230 VAC). Take note of the position of the white arrows on the plugs.
- Remove the Control Touch Pad (4) from the docking station (Fig. 7) and insert the rechargeable batteries (7) in the battery compartment on the back (Fig. 8). Make sure the batteries are inserted with the correct polarity.


**Only use the rechargeable batteries supplied or rechargeable batteries that comply with the definitions given in the Technical data. On no account should normal batteries should be used. The use of incorrect rechargeable batteries or normal batteries can cause damage to the equipment.**

- Close the battery compartment again, re-fix the Control Touch Pad to the docking station so that it clips into place, connect power supply unit (9) to the docking station (9) and plug it into the mains power (100VAC - 230 VAC). The unit is now ready for use.

**Control Touch Pad (4) with Docking Station**

- A: Touchscreen user interface
- B: „Laser in Operation“ display
- C: Laser button
- D: „Fault“ display

**Laser Unit with connected FOC**

- E: Laser Unit switch on
- F: „Fault“ display
- G: „Laser in Operation“ display

### 3. Use

3B LASER NEEDLE is a low level laser therapy for the treatment of pain-, trigger- and acupuncture points. It is intended exclusively for use in medical practices and clinics and may only be operated under permanent supervision from medical specialists.

### 4. Contraindications

The following contraindications are referred to in the laser literature:

**Absolute:**

- Irradiation of eyes
- Photosensitivity
- Tumour patients
- Open fontanelles and epiphyseal cartilaged in children
- Thyroid gland hyperactivity

**Relative:**

- Pacemaker patients (thorax)
- Epilepsy(head)
- Pregnancy (stomach and back regions)
- Endocrine organs (thorax gland, testicles, ovaries)


Therapy hindrances: deep X-ray treatment, permanent medication, chemotherapy

### 5. Operation

Switch the 3B LASER NEEDLE on at the Control Touch Pad (4) using the LASER button (C). Following the appearance of the 3B logo, the start menu opens up a few seconds later.

The first time the unit is switched on, the values 12 x 50 mW, CW and 20 minutes are pre-programmed. Otherwise, the last setting is always retained.

“Start menu” screen

1. For setting the laser power
2. For setting the operating mode
3. For setting the duration of the therapy
4. Display of the set values
5. Program saving
6. Program calling
7. Sound on/off
8. Selection of the laser diodes (red = selected)
9. Establishing wireless connection
10. Battery charge level display


“Laser in operation” screen


1. Display of the dose rate / activated laser dio
2. Display of the power level setting and the operating mode
3. Therapy duration time counts downwards
4. Activated laser diodes appear as bars, in line with the defined power setting

Select the individual parameters in the "Start menu" screen  
 POWER [mW]: Power in mW  
 PROGRAM: Operating mode  
 THERAPY TIME [min] : Therapy duration in minutes


"Power setting" screen  
 Range: 1 - 50 mW  
 Enter: Applies value to the Start menu  
 DEL: Delete key


"Operating mode setting" screen  
 Example: Standard program  
 Click on the operating mode and your selection is automatically applied to the Start menu.  
 User Defined: Program the frequencies yourself

"Operating mode setting" screen  
 Example: Nogier frequencies have been selected in the register  
 Click on the required frequency


"User Defined" screen (self-defined frequencies)  
 Range: 1.0000 to 9999.99 Hz  
 Enter: Applies value to the Start menu  
 DEL: Delete key

"Therapy time setting" screen  
 Range: 00:01 – 20:00 minutes  
 Enter: Applies value to the Start menu  
 DEL: Delete key

"Program saving" screen  
 Maximum: 40 characters / program  
 Maximum: 200 programs  
 Save  
 Home: Return to the Start menu  
 1: Defined parameters  
 2: Selected laser diodes

"Select program" screen  
 The programs are shown in alphabetical order  
 Scroll downwards with the scroll bar.  
 Home: Return to the Start menu  
 Delete: Deletes the program  
 Accept: Applies value to the Start menu


## 6. Usage

The following sequence must be followed before each laser treatment:

1. Plug the FOC selected through the Start menu into the silicone adapter. There are 2 different silicone adapters available. Always pull the FOC out of the holder individually and place them on the acupuncture point. **Only handle the FOC using the white holder.**
2. Attach the silicone adapter to the patient by applying light pressure. Fixing tape is available for the blue silicone adapter. Follow the numbering on the FOC and make sure that the numbered FOC has also been selected in the Start menu. Don't exert strong pressure or pull on the silver sleeve at the end of the FOC.

4. Once all selected FOC have been attached to the patient, switch the 3B LASER NEEDLE on at the Control Touch Pad using the LASER button.
5. Laser operation is shown on screen by "Laser in operation", at the identified LEDs and with an audible signal, every 20 seconds. The audible signal can be deactivated with the speaker button.


**Please note: No touchscreen operations are possible during laser output.**

5. 3B LASER NEEDLE switches off automatically once the defined therapy duration time has expired. Switch-off is indicated by 3 short audible signals. Treatment can be stopped at any time by pressing the LASER button. This takes you back to the Start menu.
6. Once the treatment is completed (FOC are now off), remove the individual cables from the silicone adapter starting from the left and put them back in the corresponding holders. Continue until all FOC have been dealt with.
7. Switch the 3B LASER NEEDLE off, once the treatment is completed. To do this, press and hold the LASER button for 3 seconds.
8. If there is no operator action for 30 minutes, the Control Touch Pad switches off automatically. Pressing the LASER button switches the unit back on.
9. The Laser Unit can only be switched off by unplugging the power plug.

**You may only use the 3B LASER NEEDLE if all Fiber optic cables (FOC) are connected.  
If there is no FOC on one laserdiode, the laserdiode must be protected with a cap.  
You may only exchange a FOC if the laser unit is shut off and the net plugs are disconnected.**

The Control Touch Pad can also be operated under rechargeable battery power. To do this, pull out the plug on the back of the docking station (slide the black outer ring backwards slightly).

The Control Touch Pad can also be separated from the docking station.

Please note: Once the treatment is completed, locate the Control Touch Pad back on the docking station, and click it into place. Then plug the mains power back in.

The maximum treatment time with fully charged batteries is 3 hours. The battery charging status is displayed through the battery symbol button.

**Recommendation: The Control Touch Pad should be left on the docking station at all times.  
Switch off when not in use.**

The following must be observed if you wish to achieve optimum success with the therapy.

- To keep reflection losses as small as possible, skin must be grease-free.
- Only apply creams, lotions and ampoules when the laser treatment has been finished.
- The principle of increasing application applies to all treatments. Initial treatments should not be longer than 2 minutes.

## 7. Maintenance and care

The following points must be noted and adhered to in order to protect the 3B LASER NEEDLE and in particular to prevent damage to the FOC:

**FOC should not be kinked or placed under tension. Do not tread on the FOC or pull on them. Always make sure that no knots form in the FOC and that they are not twisted. Do not apply any strong pressure to the FOC.**

**Note the minimum bend radius of 1.5 cm.**

**Only handle the FOC using the white holder and always replace them in the holders provided when not in use (Fig. 5 / Fig. 6).**

Using a soft cloth / cotton pad, clean the laser beam outlet (sleeve) using 70% ethanol following each contact with skin and when changing the patient.

- Do not expose the unit to direct sunshine.
- You can clean the housing using a gentle cleaning agent. Before cleaning disconnect net plugs.
- Humidity may not penetrate into the housing.
- Repairs must only be carried out by the manufacturer or by authorized persons. Unauthorized opening of the unit cancels the warranty. Contact the technical service if a repair is necessary.
- Only use the laser system indoors and at room temperatures

## 8. Safety information

**Please observe the instructions and warnings in the operating instructions and the warning labels on the 3B LASER NEEDLE.**

The 3B LASER NEEDLE should only be used by medical specialists who have been trained in the correct handling. Treatment must be given while under permanent supervision.

If used incorrectly, lasers can damage your eyes, especially at distances shorter than the safety distance.

**Do not look into the laser beam outlet, and do not point into the eyes of other persons even if eyes are closed.**

The safety distance when viewing directly into the beams 30 cm depending on the output power.

Make sure that no mirrors or other reflecting surfaces are located within the safety distance from the 3B LASER NEEDLE (reflection of laser light). **While treatment you always have to wear laser protection glasses (patient and therapist) for the used wavelength. Laser protection glasses have to apply to the EU-standard EN207.**

The 3B LASER NEEDLE should only be switched on once all the FOC selected in the Start menu have been attached to the patient. Once the treatment is complete, switch the 3B LASER NEEDLE off and return the FOC to their original position.

If the 3B LASER NEEDLE, the FOC, the connection plugs and/or the mains power plugs show any signs of external damage, the 3B LASER NEEDLE may no longer be operated and must be returned to the manufacturer.

Magnetic and electric fields may influence the laser system functions. Therefore do not use the LASER NEEDLE in the vicinity of equipment which generate large electromagnetic fields, e. g. diathermy units or mobile telephones.

The device and the rechargeable batteries are subject to the WEEE regulations (waste of electrical and electronic equipment) 2002/96/EG and may not be disposed off as normal household waste. If you finally wish to terminate use of the laser system, please contact the technical service .

The manufacturer is only responsible for safety and reliability if the laser system is used in accordance with the operating instructions

## Label and warning signs

### Control Touch Pad

housing upper part

laser warning sign

housing bottom part

manufacturer's label

### Laser Unit

housing front side


warning label  
outlet of laser radiation

housing rear side

warning label EN 60825-1  
labeling wavelength and laser power

housing rear side

manufacturer's label

### Identification


observe the operating  
instructions

consider disposal regulations

name of manufacturer

date of production: 2010

## 9. Accessories

W14224	3B-Dermatology-Spot
W14238	Acrylic table (60/40/90 cm)
W14239	Wall bracket
W14240	Holding arm
W14236	Fiber optic cables (FOC)
W14229	1 box with 100 single-use applicators, white
W14226	Silicone applicators, blue (pack with 5 pc)
W14227	Fixing tape (pack with 12 rolls)
W14245	Laser protection glasses - therapist (red laserlight)
W14203	Laser protection glasses - therapist (Infrared laserlight)
W14250	5 rechargeable batteries NiMH - AA Size - 2600 mAh

## 10. Status, Warning and Fault Messages

The following table lists the status, warning and fault messages and their meaning.




### 10.1 Control Touch Pad:

Element/Message	Explanation	Action
LEDs on the front panel		
Fig. 10.1.1	 Fault display The meaning is explained by the corresponding fault message on the LCD	<ul style="list-style-type: none"> <li>Check the message on the LCD</li> <li>If necessary, switch the unit off and back on</li> <li>Contact the Service team</li> </ul>
Fig. 10.1.2	 Laser output is switched on	-

Wireless interface connection status (Control Touch Pad to the Laser Unit)		
Fig. 10.1.3	 No connection to the Laser Unit	<ul style="list-style-type: none"> <li>Press the button again</li> <li>Switch the device off and unplug it from the mains</li> <li>Then plug the mains power back in and switch back on</li> <li>Contact the Service team</li> </ul>
Fig. 10.1.4	 There is a connection of the Laser Unit	

Fault messages on the LCD Confirm message with OK		
Fig. 10.1.5	 No wireless connection!	<ul style="list-style-type: none"> <li>Switch the device off and unplug it from the mains</li> <li>Then plug the mains power back in and switch back on</li> <li>Contact the Service team</li> </ul>
Fig. 10.1.6	 Low laser power detected!	<ul style="list-style-type: none"> <li>Restart laser output</li> <li>If the fault reoccurs, unplug the Laser Unit from the mains power briefly</li> <li>Then plug the mains power back in and switch back on</li> <li>Contact the Service team</li> </ul>
Fig. 10.1.7	 High laser power detected!	<ul style="list-style-type: none"> <li>Restart laser output</li> <li>If the fault reoccurs, unplug the Laser Unit from the mains power briefly</li> <li>Then plug the mains power back in and switch back on</li> <li>Contact the Service team</li> </ul>
Fig. 10.1.8	 High laser diode temperature!	<ul style="list-style-type: none"> <li>Temperature at the laser diodes &gt;60°C</li> <li>The unit switches automatically to a frequency-modulated operating mode to protect the laser diodes</li> </ul>

Rechargeable battery icons		
Fig. 10.1.9	 The battery charging level is displayed Press on the icon to display further information about the charging level	
Fig. 10.1.10	 Batteries are being charged	
Fig. 10.1.11	 Batteries not inserted or inserted with the wrong polarity	<ul style="list-style-type: none"> <li>Insert batteries into the battery compartment</li> <li>Check that the batteries are inserted the right way round</li> </ul>

Rechargeable battery status displays: Display when the battery icon is pressed			
Fig. 10.1.12	 Battery Level: 80 % State: Idle	Batteries are not being charged The battery charging level is displayed as a percentage value	
Fig. 10.1.13	 Battery Level: 80 % State: Charging ...	Batteries are being charged The battery charging level is displayed as a percentage value	
Fig. 10.1.14	 Battery Level: 0 % State: False polarity	Batteries not inserted or inserted with the wrong polarity	<ul style="list-style-type: none"> <li>Insert batteries into the battery compartment</li> <li>Check that the batteries are inserted the right way round</li> </ul>

Warning messages when saving a program Confirm message with OK			
Fig. 10.1.15	 No laser diode selected!	Saving is not possible as no laser diode has been selected	<ul style="list-style-type: none"> <li>Select at least one laser diode</li> </ul>
Fig. 10.1.16	 Incorrect name!	Saving not possible as an invalid program name has been entered (either no name or a name that already exists)	<ul style="list-style-type: none"> <li>Modify the program name</li> </ul>

Sound on / off		
Fig. 10.1.17	 Sound On	
Fig. 10.1.18	 Sound Off	

### 10.2 Laser Unit:

Element/Message	Explanation	Action
LEDs on the patient component		
Fig. 10.2.1	 Fault display The meaning is explained by the corresponding LCD on the Control Touch Pad	<ul style="list-style-type: none"> <li>Check the message on the LCD</li> <li>Unplug the Laser Unit from the mains power</li> <li>Then plug the mains power back in and switch back on</li> <li>Contact the Service team</li> </ul>
Fig. 10.2.2	 Laser Unit is switched on	<ul style="list-style-type: none"> <li>If the green LED does not light up, check whether the plug-in power supply unit is connected</li> <li>Contact the Service team</li> </ul>
Fig. 10.2.3	 Laser output active (switched on)	


## 11. Technical Data

<b>General</b>	
Product name	3B LASER NEEDLE
<b>Classification and standards</b>	
Classification:	Protection class II to EN 60601-1 Application unit BF to EN 60601-1 class IIb nach RL 93/42 EWG Laser class 3R to EN 60825-1:2007 Without moisture protection (IPX0)
<b>Laser power</b>	
Laser type:	Semiconductor laser 660/ 785/ 405 nm
Effective output CW max.	12 x 50 mW (CW Betrieb)
Wave length and choices:	<ul style="list-style-type: none"> <li>• 12 x 660 nm</li> <li>• 12 x 785 nm</li> <li>• 10 x 660 nm, 2x 785 nm</li> <li>• 11 x 660 nm, 1x 405 nm</li> </ul>
Modulation-frequencies  s. page 16	<ul style="list-style-type: none"> <li>• CW (continuous wave)</li> <li>• Multi (Multifrequency)</li> <li>• Alpha- frequency</li> <li>• Bahr frequency (1- 7)</li> <li>• Nogier- frequencies (A- G, L)</li> <li>• Reiningger frequencies 1- 12</li> <li>• Chakra frequencyies 1-8</li> <li>• A user-definable frequency can be set by the user in the range 0 to 9999 Hz</li> </ul>
Fiber optic cables	Fiber optic cables Plastic fibers Core diameter, 0.5 mm Minimum bend radius, 15 mm
Beam divergence	<ul style="list-style-type: none"> <li>• 660 nm: x-direction 7,5-12°, y- direction 15-19°</li> <li>• 785 nm: x- direction 7,5-10°, y- direction 15-19°</li> <li>• 405 nm: x- direction 6-12°, y- direction 15-21°</li> </ul>
<b>Components</b>	
Components	Control Touch Pad with Docking Station and Laser Unit
<b>RF interface</b>	
Band	ISM- Band, 2,4 GHz
Range	<= 10m
RF admission	to EN 300440-1/2

<b>Control Touch Pad + Laser Unit</b>	
Power supply	Power supply unit Optional with rechargeable batteries
	Primary: 100-240 VAC/ 50-60 Hz/ 400 mA
	Secondary: 12 V/1,25 A/ 12 W, short-circuit-proof
	Manufacturer: Friwo, Typ FW 7555M/12
	Application unit to EN 60601-1
Fuses	2 A , power supply unit short-circuit-proof
Rechargeable batteries	5 pc, AA- Size, NiMH, 1,2 V Manufacturer f.e. Panasonic , HHR-3XRE, 2600 mAh The service life of a full charge > 4 h Charging duration: approx 5h
<b>Permissible temperature range</b>	
Permissible temperature range	10 - 30 Grad C
Storage temperature	0 - 50 Grad C
Permissible relative humidity	30 – 75 %
<b>Manufacturer</b>	
Manufacturer	livetec GmbH, Marie- Curie- Str. 8, 79539 Lörrach, Germany
<b>Conformity</b>	
Conformity	CE 0366
<b>WEEE-Reg.- Nr.</b>	
WEEE-Reg.- Nr.	FE 59335168
<b>Exclusive sales</b>	
	3B Scientific GmbH, Rudolfweg 8, 21031 Hamburg, Germany Phone: +49-(0)-40-73966-222 Fax: +49-(0)-40-73966-100 e-mail: stefan.baudis@3bscientific.com Internet : www.3b-laser.de

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Hersteller: **Livetec GmbH, Marie-Curie-Str. 8, 79539 Lörrach - Germany**
Lieferant: 3B Scientific GmbH