








E-SCOPE®

HEARING IMPAIRED



-  Used with hearing aid or hearing impairment
-  One-handed control over volume adjustments
-  Specialist adult diaphragm
-  Automatic shut off after 1.5 to 2 minutes
-  Audio output jack
-  Battery operated (single AAA battery)
-  Supplied with accessory pack



The **E-Scope Hearing Impaired Model** includes an assortment of auditory accessories, allowing users with hearing loss the ability to choose the appropriate listening device for their specific needs. This model is identical in functionality to the Clinical Model — amplifying sounds up to 30 times louder than an acoustic scope. Separate headphones and accessories are offered to ensure we find the right solution for any user with hearing loss.

*E-Scope Catalog # 718-7700 (Clinical model)
718-7710 (Hearing impaired model)
Technical Support: Unlimited
Warranty: One year – parts and labor*

Cardionics
™ The Heart of Auscultation 

Phone: 281-488-5901, 800-364-5901 (USA)
Fax: 281-488-3195
Email: info@cardionics.com
Web: www.cardionics.com
Address: 910 Bay Star Blvd, Webster, TX 77598 USA