



## Eye-Level Digital Physician Scales

**500 lb x 0.2 lb / 225 kg x 0.1 kg**

These 59" (150 cm) high digital scales features the eye-level "at-a-glance" convenience of a battery-powered 0.5"/12.7 mm high-contrast LCD display. With the push of a button, readings appear in pounds OR kilograms and the weight display also calculates **Body Mass Index**. The 10.5" W x 14.5" D (27 cm W x 37 cm D) platform cover is easily removed for cleaning. Handpost and height rod combinations are also available. Please see the back page for different model numbers.

Detecto's eye-level digital physician scales feature the powerful 758C digital weight indicator, which has battery or AC power, patient ID storage, serial connectivity for output to a printer or PC, and a host of other useful functions.

**YOUR CHOICE: MECHANICAL OR DIGITAL HEIGHT RODS**

# Detecto Digital Physician Scale Features



11:23, 10/15/05  
 ID: 12548796  
 BMI: 28.4  
 Height: 5' 11.0"  
 Weight: 204 lb G



**With the Enhanced 758C  
 Digital Weight Display**

## SERIAL CONNECTIVITY PROVIDES:

**OUTPUT TO PRINTER** print weight, BMI, time, date, and height

**SEND DATA TO PC** excellent for efficient electronic medical record keeping

### DETECTO'S EYE-LEVEL PHYSICIAN SCALE FEATURES THE 758C WEIGHT DISPLAY

Pushbutton and keypad tare, automatic display lock, and selectable units make the 758C the perfect choice for weight loss clinics and long-term care facilities.

### BODY MASS INDEX

Body Mass Index (BMI) calculations have been shown to be an accurate estimate of health risk when used as part of a comprehensive health assessment. The 758C quickly calculates BMI by simply entering the patient's height while weighing.



### SERIAL CONNECTIVITY

The 758C is the new standard for health care scale weight indicators. Serial connectivity (illustrated above) provides output to a PC or printer for efficient record keeping, useful for patient records in physicians' offices, weight loss clinics, or athletics. The new serial port allows output of the 758C's data directly to a printer by simply pressing the PRINT/ENTER key. The paper copy can be included in a patient's files.



### PATIENT ID

Easily enter up to an 11-digit patient ID and print it along with the weight, time and date from a label printer or send the data to a PC.

MODEL	CAPACITY	DESCRIPTION
6437	500 lb x 0.2 lb / 225 kg x 0.1 kg	No Height Rod
6437DHR	Same	Digital Height Rod
6439	Same	Mechanical Height Rod
6449	Same	Handpost / Mech. Ht. Rod

**Model 6439/6449 Height Rod:** 30" - 78" and 76 cm - 200 cm

*Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.*



The 758C is **battery powered** for use with portable scales, so you can take the scale to the patient. Optional AC adapter available.



**Sold By:**

**DETECTO**  
 A Division of Cardinal Scale Manufacturing Co.  
 203 E. Daugherty, Webb City, MO 64870 USA  
 Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001  
 www.Detecto.com