

PARALLEL BARS



AM-700 Parallel Bars - Electric/Platform Mount 10' ⁷ also available AM-702

- Height: 26" - 48"
- Width: 15" - 29"

Power height-adjustable parallel bars. Digital Display on conveniently located user-friendly control panel shows exact bar height. Memory buttons allow storage and recall of four separate bar height settings. Microprocessor control unit ensures synchronous operation of all four columns. Heavy-gauge 1 1/2" diameter (3.8 cm) stainless steel handrails on three-section slim line lifting columns that provide maximum lateral stability, even in the highest bar position. Height monitoring functions ensure reliable operation and safe switch-off in case of a fault.

- Slip resistant laminated platform
- Platform 10'(3m)L x 43"W
- 29"W between columns
- Stainless steel handrails 1 1/2" diameter (3.8cm)
- Height adjusts electrically from 26" to 48"
- Width adjusts manually from 15" to 29"
- Weight capacity on hand rails: 400 lbs.